

Ensemble CMC™

Anatomy Saving By Design

Thumb arthritis?

Ask about the minimally-invasive solution



Ensemble CMC™

***Has your doctor
recommended joint
replacement surgery
for treatment of mild to
moderate osteoarthritis
of the thumb?***

A minimally invasive option is now available from Ensemble Orthopedics that may help reduce pain and improve joint function.

The Ensemble CMC™ is an anatomically shaped, stemless device that is inserted through a small incision into the space between the damaged articular surfaces in the CMC joint of the thumb. The minimally invasive surgical procedure does not disturb the stabilizing soft tissues of the joint, allowing potential for a shorter recovery time.

The Ensemble CMC™ interpositional implant is made from pyrocarbon – a biocompatible, low-friction, and durable material with a long, proven history of use in the manufacture of artificial heart valves and orthopedic devices.

For more information, scan this code with your smartphone:



www.ensembleortho.com
info@ensembleortho.com
+1 (512) 638-3598

MK-72-101-02 rev C